

130mm

Singclean®



# COVID-19 IgG/IgM Test kit (Colloidal Gold Method)

**INTENDED USE**

COVID-19 IgG/IgM Test kit (Colloidal Gold Method) is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus in human whole blood, serum or plasma. This test provides only a preliminary test result. Therefore, any reactive specimen with the COVID-19 IgG/IgM Test kit (Colloidal Gold Method) must be confirmed with alternative testing method(s) and clinical findings.

**PACK FORMATS**

1 test/box ,20 test/box,50 test/box , 100 test/box

**INTRODUCTION**

Coronaviruses are enveloped RNA viruses that are distributed broadly among humans, other mammals, and birds and that cause respiratory, enteric, hepatic, and neurologic diseases. Seven coronavirus species are known to cause human disease. Four viruses - 229E, OC43, NL63, and HKU1 - are prevalent and typically cause common cold symptoms in immunocompetent individuals. The three other strains - severe acute respiratory syndrome coronavirus (SARS-CoV), Middle East respiratory syndrome coronavirus (MERS-CoV) and 2019 Novel Coronavirus (COVID-19) - are zoonotic in origin and have been linked to sometimes fatal illness. IgG and IgM antibodies to 2019 Novel Coronavirus can be detected with 2-3 weeks after exposure. IgG remains positive, but the antibody level drops overtime.

**PRINCIPLE**

The COVID-19 IgG/IgM Test kit (Colloidal Gold Method) uses the principle of colloidal gold immunochromatography. The test uses anti-human IgM antibody (test line IgM), anti-human IgG (test line IgG) and goat anti-mouse IgG (control line C) immobilised on a nitrocellulose strip. The burgundy colored conjugate pad contains colloidal gold conjugated to recombinant COVID-19 antigens conjugated with colloidal gold (COVID-19 conjugates). When a specimen followed by assay buffer is added to the sample well, IgM &/or IgG antibodies if present, will bind to COVID-19 conjugates making antigen antibodies complex. This complex migrates through nitrocellulose membrane by capillary action. When the complex meets the line of the corresponding immobilized antibody (anti-human IgM &/or anti-human IgG) the complex is trapped forming a burgundy colored band which confirm a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result.

Regardless of the presence of COVID-19 IgG/IgM antibody in the test sample, a colored band will appear on the quality control line (line C). The colored band appearing on the quality control line (line C) is used to determine whether there are enough samples and whether the chromatographic process meets the normal standard, and also serves as the internal control standard of the reagent.

**MATERIALS SUPPLIED**

- Sealed pouches each containing a test cassette, a desiccant dropper
- Lancets (for fingerstick whole blood only)
- sterilizing tablet (for fingerstick whole blood only)
- Buffer
- Package insert

**MATERIAL REQUIRED BUT NOT PROVIDED**

1. Specimen collection containers
2. Centrifuge (for plasma only)
3. Timer

**STORAGE AND STABILITY**

The kit can be stored at room temperature or refrigerated (4-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

**WARNINGS AND PRECAUTIONS**

1. For In Vitro diagnostic use only. Do not use after expiration date.
2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
3. Do not use it if the tube/pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, i.e. electric fan or strong air-conditioning.
9. COVID-19 IgG/IgM Test Kit (Colloidal Gold Method) (hereinafter referred as Product) used to detect antibodies produced by human body to against COVID-19 after infected. COVID-19 positive patients may be tested negative by this Product at that period, as human body has not produced or produced very few antibodies at that period, which can't be detected by the Product.
10. Product can be only used as auxiliary diagnosis for patients infected with COVID-19 and should not be the final diagnosis of COVID-19 infection. It must be combined with other medical tests for comprehensive diagnosis.
11. The Product is designed for use by healthcare professionals. It should be operated by healthcare professionals. Individuals are not allowed to use.

**SPECIMEN COLLECTION**

1. COVID-19 IgG/IgM Test kit (Colloidal Gold Method) can be performed using either whole blood, serum or plasma.
2. Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
3. Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
4. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
5. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

**TEST PROCEDURE**

Allow test cassette, specimen, buffer and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test device on a clean and level surface.

**For Serum or Plasma Specimens:**

With the provided dropper, draw serum/plasma specimen and then add 1 drop serum/plasma specimen into the sample well(S). Then add 2 drops of buffer to the sample well(S) immediately. Avoid air bubbles.

**For Whole Blood Specimen:**

With the provided dropper and transfer 1 drop of whole blood to the sample well(S) of the test device, then add 2 drops of buffer to the sample well(S) immediately. Avoid air bubbles.

3. Wait for the colored line(s) to appear. The result should be read after 10 minutes. Do not interpret the result after 20 minutes.

**OPERATION INSTRUCTIONS**

correct use, convenient and quick



1. Constant Temperature and Static position;



2. Massage finger pulp;



3. Disinfect with alcohol tablet;




4. Blood sampling and acupuncture into the abdomen of fingers;



5. Wipe off the first drop of blood with sterile dry cotton ball;



6. Use disposable micro pipette to suck blood;



**INTERPRETATION OF RESULTS**

**NEGATIVE:** If only the C band is present, the absence of any burgundy color in the both T bands (IgG and IgM) indicates that no anti-COVID-19 antibodies are detected in the specimen. The result is negative.

**IgM POSITIVE:**

In addition to the presence of C band, if only IgM band is developed, the test indicates for the presence of IgM anti-COVID-19 in the specimen. The result is IgM anti-COVID-19 positive.

**IgG POSITIVE:**

In addition to the presence of C band, if only IgG band is developed, the test indicates for the presence of IgG anti-COVID-19 in the specimen. The result is IgG anti-COVID-19 positive.

**IgG and IgM POSITIVE:**

In addition to the presence of C band, both IgG and IgM bands are developed, the test indicates for the presence of both IgG and IgM anti-COVID-19 in the specimen. The result is IgG and IgM anti-COVID-19 positive.

**INVALID:**

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**QUALITY CONTROL**

A procedural control is included in the test. A red line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

**LIMITATIONS**

1. Use fresh samples whenever possible. Frozen and thawed samples (especially repeatedly) contain particles that can block the membrane. This slows the flow of reagents and can lead to high background color, making the interpretation of results difficult.
2. Optimal assay performance requires strict adherence to the assay procedure described in this insert sheet. Deviations may lead to aberrant results.
3. A negative result for an individual subject indicates absence of detectable anti-COVID-19 antibodies. However, a negative test result does not preclude the possibility of exposure to or infection with COVID-19.
4. A negative result can occur if the quantity of the anti-COVID-19 antibodies present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
5. Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
6. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

**PERFORMANCE CHARACTERISTICS**

1. Clinical Performance for IgM Test  
The samples from susceptible subjects were tested by The COVID-19 IgG/IgM Test kit (Colloidal Gold Method) and by a commercial IgM EIA kit.  
Relative Sensitivity: 95.7%, Relative Specificity: 97.3%, Overall Agreement: 96.8%
2. Clinical Performance for IgG Test  
The samples from susceptible subjects were tested by The COVID-19 IgG/IgM Test kit (Colloidal Gold Method) and by a commercial IgG EIA kit.  
Relative Sensitivity: 91.8%, Relative Specificity: 96.4%, Overall Agreement: 95.0%

**REFERENCE**

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85-164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology. 6th ed. Lippincott Williams & Wilkins, 2013: 825-58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490-502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181-192.

**SYMBOLS USED ON PACKAGING NOTICE**

	Authorized Representative		Store between 4-30°C		For in vitro diagnostic use only
	Do not reuse		Lot Number		Don't use if package is damaged
	Consult instructions for use		Use by		

331mm